



**Ride Designs®**  
 a branch of Aspen Seating, LLC  
 toll-free 866.781.1633  
 phone 303.781.1633  
 fax 303.781.1722  
 www.ridedesigns.com

## Ride Forward™ Cushion Order Form

<b>Supplier Name &amp; Address</b>	<b>Ship To</b>	<b>Details</b>
_____	_____	Account # _____
_____	_____	P.O. # _____
_____	_____	Purchaser _____
Phone _____	Phone _____	Mark for _____
Fax _____	Fax _____	Date _____

### Ride Forward Cushion (E2607, skin protection and positioning wheelchair seat cushion)

Forward Cushion Size (width x depth)	Part Number Complete	Mfr. Sugg. Retail Price* Complete	Fits Wheelchair Widths	Fits Wheelchair Depths	Quantity	Total
12 x 12 <b>P</b>	FC-1212	\$425.00	11/12"	12/13"		
14 x 14	FC-1414	\$425.00	13/14"	14/15"		
14 x 16	FC-1416	\$425.00	13/14"	16/17"		
16 x 16	FC-1616	\$425.00	15/16"	16/17"		
16 x 18	FC-1618	\$425.00	15/16"	18/19"		
16 x 20	FC-1620	\$425.00	15/16"	20/21"		
18 x 16	FC-1816	\$425.00	17/18"	16/17"		
18 x 18	FC-1818	\$425.00	17/18"	18/19"		
18 x 20	FC-1820	\$425.00	17/18"	20/21"		
20 x 16	FC-2016	\$425.00	19/20"	16/17"		
20 x 18	FC-2018	\$425.00	19/20"	18/19"		
20 x 20	FC-2020	\$425.00	19/20"	20/21"		

Subtotal Ride Forward Cushions \_\_\_\_\_

### Additional Forward Items

Item	Part Number	Mfr. Sugg. Retail Price*	Quantity	Total
Outer Cover Only -- specify Cushion width x length --	FCE-_____	\$125.00	_____	
Ride CAM® Wedge Kit* -- specify Cushion width --	FCC-WK _____	\$ 50.00	_____	

Subtotal Additional Forward Items \_\_\_\_\_

**TOTAL from all categories** \_\_\_\_\_

\* NOTE: Each kit contains two (2) Ride CAM wedges for 14" to 20" width cushions only.

**P** = Contoured and sized for pediatric application.

### Return Policy

Ride Forward Cushions may be returned if they are in their original, sealed, Ride plastic bag. All returns are subject to a 10% restocking fee.

\*All prices are in U.S. dollars.

### Have you considered a back support?



### The Corbac™ Adjustable Contour Back Support:

- lightweight
- breathable
- easy adjustment
- no hardware
- folds with wheelchair

E2611, General use wheelchair back cushion

### Ride CAM® Wedge Kit



Use Ride CAM Wedges **symmetrically** (shown above) to promote more upright posture and decrease ischial pressure or **asymmetrically** to help level a flexible pelvic obliquity.

Each kit contains two (2) Ride CAM wedges for 14" to 20" width cushions only.