

New!

Ride[®] Custom Cushion 2

Announcing the marriage of a new Ride Custom Cushion + advanced scanning/ordering technology!

Fast, clean, intuitive and accurate process.

Rapidly scan, then send shape and order information through the secure, encrypted Rideworks™ iOS app.

Touch “send” and we’ll do the rest... providing a custom cushion made specifically to match each individual’s unique shape, posture, and skin care needs.

Ride Custom Cushion 2: All the skin and postural care of the original Custom Cushion — now in a more comfortable, lighter and simpler design.

It's all about contour...



Ride Custom Cushion 2

▲ Spacer fabric in cover.



Helps protect skin and enhance posture

Ride's patented and evidence-based mechanism of support helps to enhance sitting stability and performance, while reducing the dangerous forces of pressure and shear from at-risk bony prominences.

Lightweight

About 1.8 pounds*: A fraction of the weight of other cushions.

*Weight of an average 16" by 16" size Ride Custom Cushion 2.

Breathable

Contours and spacer fabric cover promotes airflow to help keep the sitter dry and comfortable in virtually any climate.

Simple and adjustable

The Ride CAM[®] wedges can be used to adjust the contours to accommodate weight change, tissue atrophy, and postural change.

Growable

Cushion can be grown in width and length at a fraction of the cost of a new cushion.

RIDEWORKS™

NEW Simulation, Scanning, & Ordering Process



No more plaster or impression foam simulators!

Ride Custom Cushion 2 and Ride Custom Back simulation can now be completed with a dual bead-bag and vacuum system. The patent-pending Shape Capture Base, coupled with the seat bead-bag, helps to stabilize the pelvis while the bead bag is malleable, and enhances shape accuracy. It also creates the necessary reference for the seat contour orientation. No need for laser levels, putty or roping.

Fast, accurate and safe

Capture shapes, draw trim-lines directly on the bead bags and then let the new RideWorks scanner and app guide you through the rest. A fresh outer simulator bag for each simulation helps prevent cross-contamination and infection.

RideWorks

In minutes, Ride's exclusive app organizes the process of client photos, scanning, and order form completion — all of which are then transmitted to Ride Designs via the internet with secure encryption.

Use with Confidence

The Ride Custom Cushion 2 comes with a 90 day fit and function guarantee and an 18 month limited warranty.


The Cushion has a maximum client weight of 250 pounds for cushions 20" width or narrower; 300 pounds for 21" or wider.

Make it a System!




Pair a Ride Custom Back with the Ride Custom Cushion 2 for unparalleled sitting performance

Better Outcomes



The Old: Immersion/Envelopment

Pressure map of a person on a properly-inflated air style cushion. Note peak pressure circled at right ischial tuberosity.



The New: Ride Designs

Same person on a Custom Cushion by Ride Designs. Note off-loading of both ischial tuberosities (A), coccyx/sacrum (B), and decreased pressure at greater trochanters (C).

The Ride Custom Back

Lightweight strength and accuracy Thin profile has little or no impact on wheelchair seat depth, width, or performance, making the Custom Back appropriate even for active users.

Adjusts easily The infinitely adjustable FlexLoc® mounting system allows adjustment in every axis of movement, ensuring optimal balance and postural control.

Breathable The Brock Composite liner material and spacer fabric cover help to keep the sitter dry and comfortable in virtually any climate.

Easy to clean Removable and ventilated inner shell eases cleaning and maintenance.

Growable Like the Ride Custom Cushion 2, the Back can be grown at a fraction of the cost of a new custom system.



Ride Designs®
a branch of Aspen Seating, LLC
8100 SouthPark Way, C400
Littleton, Colorado 80120 USA

phone: 303.781.1633
toll-free: 866.781.1633
www.ridedesigns.com